Ultimaster™ Tansei™

Ultimaster Tansei Sirolimus eluting coronary stent system
Ultimaster Tansei, système de stent coronaire à libération de sirolimus
Sirolimus-beschichtetes Koronar-Stentsystem Ultimaster Tansei
Sistema de Stent Coronario Liberador de Sirolimús Ultimaster Tansei
Sistema de Stent Coronário com Eluição de Sirolimus Ultimaster Tansei
Sistema di Stent Coronarico Ultimaster Tansei a Rilascio di Farmaco Sirolimus
Ultimaster Tansei Sirolimus-eluerend Coronair Stentsysteme
Läkemedelseluerande Koronarstentsystem Ultimaster Tansei Sirolimus
Ultimaster Tansei Sirolimus lægemiddeleluerende Koronar-stentsystem
Ultimaster Tansei Sirolimusfrigivende Koronar-stentsystem
Ultimaster Tansei Sirolimus-Sepelvaltimostenttijärjestelmä
Συστημα Στεφανιαιασ Ενδοπροθεσησ με Εκλυση Φαρμακου Sirolimus Ultimaster Tansei
Сиролимус-выделяющий Коронарный Стент Ultimaster Tansei
System Stentu Wieńcowego Uwalniającego Lek Sirolimus Ultimaster Tansei
Ultimaster Tansei Sirolimus-kibocsátó Koronária Sztent Rendszer
Koronární stentový systém Ultimaster Tansei uvolňující léčivo Sirolimus
Systém Koronárneho Stentu Ultimaster Tansei uvoľňujúci Sirolimus
Ultimaster Tansei Sirolimus Salan Koroner Stent Sistemi
Ultimaster Tansei Sirolimus Ravimkaetud Pärarteri Stendisüsteem
Ultimaster Tansei Sirolimus Medikamēntu Izdaloša Koronārā Stenta Sistēma
„Ultimaster Tansei Sirolimus“ Vaistus Išķiriančio Koronarinio Stento Sistema
Sistem Koronarne Opornice Ultimaster Tansei Sirolimus, ki sprošča zdravilo
Ultimaster Tansei Sirolimus Koronarni Stent Sistem Koji Oslobađa Lek
Sistemul Stentului Coronarian Cu Eliberare de Sirolimus Ultimaster Tansei
Коронарни стентове Ultimaster Tansei с отделяне на сирулимус
Стент-система коронарна сіролімус-елютуюча Ultimaster Tansei
Sustav koronarnog stenta za izlučivanje sirolimusa Ultimaster Tansei


DE34M302-02

REF

LOT

STERILE R


NP

Catalogue number Référence du catalogue Bestellnummer Número de catálogo Referência Codice prodotto Catalogusnr. Artikelnummer Varenummer Artikelnummer Tuotekoodi Κωδικός Προϊόντος Каталожный номер Kod produktu Termékkód Kód výrobku Kód výrobku Ürün kodu Kood Koda Nr. Katalogo Nr. Koda izdelka Šifra proizvoda Codul produsului Каталожен номер Номер за каталогом Kataloški broj	Batch code Numéro de lot Chargenbezeichnung Número de lote Referência Numero di lotto Lotnummer Batchkod Batchnummer Batch nummer Eränumero Αριθμός παρτίδας Номер партии Uzýçý przed datą Tételszám Číslo výroby serie Identifikačné číslo Lot numarası Partii nr. Seriñas Nr. Partijos kodas Številka serije Serijska oznaka Numar lotului Партиден номер Нод партиі Šifra serije	Use by date A utiliser jusque Verwendbar bis Caducidad Data limite de uso Data di scadenza Gebruiken voor Användes före Udløbsdato Brukes for-dato Eräpäivä Ημερομηνία λήξης Годен до Uzýçý przed datą Felhasználható Použití do Použití do Son kullannin tarihi Säilivusaeg Deriguma terminš Sunaudoti iki Uporabno do Upotrebljivo do Data de expirare Срок на годност Використати до Upotrijebiti do	Do not use if package is damaged Ne pas utiliser si l'emballage est endommagé Inhalt beschädigter Packung nicht verwenden No usar si el paquete está dañado Não usar se a embalagem estiver danificada Non utilizzare se la confezione risulta danneggiata Niet gebruiken wanneer de verpakking beschadigd is Använd inte om förpackningen är skadad Må ikke anvendes, hvis pakningen er beskadiget Må ikke brukes hvis forpakningen er ødelagt Älä käyttää jos pakkaus on vahingoittunut Μην το χρησιμοποιείτε εάν η συσκευασία είναι καταστραφεί Не использовать, если упаковка повреждена Nie stosować w przypadku uszkodzenia opakowania Ne használja, ha a csomagolás sérült! Непови́жайте, бы-ло-ли ба-лені пошкоже-но Непови́жайте, ак је обал пошкоде-но Paket használása kullannymiz Ärge kasutage, kui pakend on kahjustatud Nelleto, ja steriilais iepakungas ir bojāts Jei pakotėje pažeista, nenaudoti Ne uporabljajte ga, če je ovojnina poškodovana. Ne upotrebjavati ukoliko je pakovanje oštećeno A nu se utiliza în cazul în care ambalajul este deteriorat Не използвайте, ако опаковка е повредена не використовувати, якщо пакування пошкоджено Ne upotrebjavati ako je pakiranje oštećeno	Sterilized using irradiation Stérilisé par irradiation Sterilisation durch Bestrahlung Esterilizado utilizando irradiación Esterilizado com irradiação Sterilizzato con radiazioni ionizzanti Gesteriliseerd door bestraling Strålningssteriliserad Strålesteriliseret Strålesteriliseret Steriloitu säteilyttämällä Αποστειρωμένο με ακτινοβολία Стерилизовано облучением Sterilizowane przy użyciu promieniowania Besugárzással sterilizálva Sterilizovano ozářováním Sterilizované žiarením Irradysayon ile sterilize edimıştır Steriliseeritud radiatsiooniga Sterilizēts ar irradiāciju Sterilizovana jonizujućom spindulovom Sterilizirano z elektronskim žarkom Sterilisanu radijacijom Sterilizare prin radiație Стерилизован с радиация Простерилизовано за допомогою випромінювання Sterilizirano zračenjem	Do not resterilize Ne pas re-stériliser Nicht resterilisieren No reesterilizar Não reesterilizar Non risterilizzare Niet opnieuw steriliseren Får ej omsteriliseras Må ikke reesteriliseres Må ikke reesteriliseres Ei saa uudelleen steriloida Μην το επαναποστειρωθεί Не стерилизовать повторно Nie sterylizować ponownie Resterilizációja nem megengedett! Opakované nesterilizujte Nesterilizovať Sterilizujte etmeyin Ärge steriliseerige uuesti Nesterilizēt atkārtoti Nesterilizuoti Ni namerjeto ponovni sterilizaciji Ne reesterilisiati Nu resterilizați Да не се стерилизира повторно Повторно не стерилизувати Ne ponovo sterilizirati	Do not reuse Ne pas réutiliser Nicht zur Wiederverwendung No reutilizar Não reutilizar Monouso Niet opnieuw gebruiken Ingen återanvändning Må ikke genbruges Må ikke gjenbrukes Kertakäyttöinen Μόνο για μια χρήση Не использовать повторно Nie używać powtórnie Csak egyszeri használatra Nelze použít opakovaně Pro jednorázové použitie Tek kullanılmık Ainult ühekordseks kasutuseks Tikai vienkreizējai lietošanai Pakartotina nenaudoti Samo za jednokratno uporabo Samo za jednokratnu upotrebu Doar pentru unica folosinta Саме за еднократно употреба Використувати повторно заборонено Ne upotrebjavati ponovno	Temperature limitation for storage Limites de température pour le stockage Temperaturbegrenzung für die Lagerung Limitación de la temperatura de almacenamiento Limitação de temperatura para armazenamento Limitazione della temperatura per la conservazione Temperatuurimieten tijdens bewaren Temperaturbegrensning för förvaring Temperaturgrænser under opbevaring Temperaturgrænser for lagring Lämpötilarajotus Περιορισμός θερμοκρασίας για την αποθήκευση Предельные температуры хранения Ograniczenia temperatury do przechowywania Tárolási hőmérséklet korlátozás Teplotní pásmo skladování Teplotné obmedzenie pre skladovanie Saklama için sıcaklık kısıtlaması Temperatuuripiirang säilitamisel Temperatūras ierobežojumi uzglabāšanai Laikymo temperatūros ribos Temperatūras ierobežojumi za shranjevanje Dozvoljena temperatura pri likom skladištenja Limita de temperatură Ограничение на температурата за съхранение Допустима температура зберігання Temperaturno ograničenje za pohranu	Nominal Pressure Pression nominale Nenndruck Pressão nominal Pressão nominal Pressione nominale Nominale druk Nominell tryck Nominell tryk Nominell trykk Nominiaalipaine Ονομαστική Πίεση Номинальное давление Ciśnienie znamionowe Névleges nyomás Nominální tlak Menoviti tlak Nominal Basıncı Nominiaalrõhk Nominālais spiediens Minimalus slėgis Nominální tlak Nominalni pritisk Presiune nominală Ограничение на давление Номинальний тиск Nazivni tlak
---	---	---	---	---	--	---	---	--

RBP

P

Stent I.D./ Ø <L>

Contents

Manufacturer

Manufacturing site

MR conditional

Rated Burst Pressure Pression nominale de rupture Zerreißdruck Presión de rotura Pressão nominal de ruptura Pressione massima di sicurezza Bepaalde scheurdruk Nominell bristningstryck Nominel øvre trykgrænse Nominell sprengtrykk Nimellinen rikkoutuspaine Ονομαστική Πίεση Ρήξης Расчетное давление разрыва Znamionowe ciśnienie rozrywające Becsült szakítónyomás Maximální tlak Menoviti tlak prasknutia Anma Patlama Basıncı Arvestuslik purunemisrõhk Aprēķinātais sprādziena spiediens Minimalus išplėtimo slėgis Naznačen delovni tlak Pritisak Presiunea nominală de explozie Номинално налягане на разкъсване Номинальний тиск розриву Nazivni tlak pucanja	Pressure Pression Druck Presión Pressão Pressione Druk Tryck Tryk Paine Πίεση Давление Ciśnienie Nyomás Tlak Tlak Basiņç Rõhk Spiediens Slėgis Tlak Pritisak Presiune Налягане Тиск Tlak	Guiding Catheter minimum inner diameter Diamètre interne minimum du cathéter guide Mind. Innendurchmesser Führungskatheter Diámetro interior mínimo del catéter guía Diámetro interno mínimo cateter guia Minimum binnendiameter van de guiding catheter Ledarkateter minsta innerdiameter Styrekateter mindste indvendige diameter Føringskateter minimum innvendig diameter Ohjainkatetrin pienin sisähalkaisija Ελαχιστή εσωτερική διάμετρος Οδηγού Καθετήρα Минимальный внутренний диаметр проводникового катетера Średnica wewnętrzna stentu Vezető katéter minimum belső átmérő Tlak Minimální vnitřní průměr zaváděcího katétru Vodiaci katéter - minimálny vnútorný priemer Kilavuz Kateter minimum iç çapı Juhtekatetri minimaalne sisediameter Vaditāj katetra minimālais iekšējais diametrs Kreipiamojo kateterio minimalus vidinis skersmuo Najmanjši notranji premer vodilnega katetra Minimalni unutrašnji prečnik uvodnog katetera Diametru interior minim al cateterului de ghidare Минимален вътрешен диаметър на въвеждащия катетър Мінімальний внутрішній діаметр направляючого катетера Kateter za uvođenje, minimalni unutarnji promjer	Stent inner diameter Diamètre interne du stent Innendurchmesser des Stents Diámetro interior del Stent Diámetro interno do Stent Binnendiameter van de stent Stentens Innerdiameter Innvendig diameter Føringskateter på stenten Stentín Sisähalkaisija Εσωτερική διάμετρος της ενδοπρόθεσης Внутренний диаметр стента Średnica wewnętrzna stentu Szent belső átmérő Vnitřní průměr stentu Vnútorný priemer stentu Stent iç çapı Stenti sisediameter Stenta iekšējais diametrs Stento vidinis skersmuo Notranji premer opronice Unutrašnji prečnik stenta Diametru interior al stentului Вътрешен диаметър на стента Внутрішній діаметр стенту Unutarnji promjer stenta	Stent length Longueur du stent Länge des Stents Longitud del Stent Comprimento do Stent Lunghezza Stent Lengte van de stent Stentens längd Stentens længde Lengden på stenten Stentín pituus Μήκος της ενδοπρόθεσης Длина стента Długość stentu Szent hossz Délka stentu Dĺžka stentu Stent uzunluğu Stenti pikkus Stenta garums Stento ilgis Dołżina opronice Dužina stenta Lungimea stentului Дължина на стента Довжина стенту Dužina stenta	Consult instructions for use Lire le mode d'emploi Gebrauchsanweisung beachten Consulte las instrucciones de uso Consultar instruções de uso Consultare le istruzioni per l'uso De gebruiksaanwijzing raadplegen Läs bruksanvisningen Læs bruksanvisningen Les bruksanvisning for bruk Lue käyttöohjeet Διαβάστε τις οδηγίες χρήσης Прочитайте инструкции по применению Sprawdź w instrukcji stosowania Olvassa el a használati útmutatót! Postupujte podle pokynů Prečítajte si návod na použitie Kullanma talimatına başvurunuz Tutvu kasutusjuhendiga. Pirms lietošanas izlasiet instrukciju Skaitykite naudojimo instrukciją Glej navodila za uporabo Prečitate instrukcije za upotrebu A se consulta instrucțiunile de folosință Прочетеете указанията за употреба Ознайомте із інструкціями для застосування Pogledati upute za upotrebu	Contents Contenu Inhalt Contenido Conteúdo Contenuto Inhoud Innehåll Indhold Innhöld Sisältö Perleχόμενα Содержание Zawartość Tartalom Obsah Obsah İçindekiler Sisu Satus Sudėtis Vsebina Sadržaj Съдържание Зміст Sadržaj	Manufacturer Fabricant Hersteller Fabricante Fabricante Fabrikant Tilverkare Produceret af Producent Valmistaja Παράσκευαστής Производитель Wytwórca Gyártó Vyrobeno Výrobca Üretici firma Tootja Räätöläjs Gamintojas Izdevalolec Sadržaj Fabricat de Произведено от Виробник Proizvođač	Manufacturing site Site de fabrication Herstellungsort Lugar de fabricación Site do Fabricante Sito di produzione Productielocatie Tilverkningsplats Produktionssted Produksjonssted Valmistuspaikka Τόπος Κατασκευής Произведено Miejsce produkcji Gyártóhely misto výroby Miesto výroby Üretici firma Tootmiskoht Räätöläjs Ražotnes vieta Proizvodno mesto Mesto proizvodnje Local de Productie масто на производство Виробнича дільниця Mjesto proizvodnje	MR conditional Compatible avec IIRM sous certaines conditions MR-tauglich RM condicional Condicionato a RM Compatibilità alla RM con riserva MR-conditioneel MR-vilkorlig MR-betinget MR conditional MR-ehdollinen Συμβατή με MR Условно совместим с MPT Warunkowa zgodność z MRI MR-feltételes Podmíněné bezpečný při vyšetření MRI Podmínečně bezpečný v prostředí MR MR şartı RI tingimustele vastav MR nosacijums Magnetinio rezonanso sąlygos Pogojna uporaba pri slikanju MR MR-ehdollinen Compatibilitate RM conditionată MR условен Умовно сумісний з MPT Uvjeto podnosi MR
--	---	--	--	--	--	---	--	--	---

PRODUCT DESCRIPTION

Ultimaster Tansei Sirolimus Eluting Coronary Stent System with Rapid Exchange Balloon Delivery System consists of a balloon expandable intra-coronary L605 cobalt chromium (CoCr) stent with abluminal drug eluting coating, that consists of a blend of Sirolimus and poly (D,L-lactide-co-caprolactone), pre-mounted onto a high pressure, semi-compliant balloon delivery catheter.

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
DEVICE COMPONENT DESCRIPTION

STENT PLATFORM

The platform for Ultimaster Tansei Sirolimus Eluting Coronary Stent System, designed based on the Kaname® CoCr Coronary Stent, offers a flexible stent with thin struts, excellent deliverability and optimized side branch access.

STENT COATING LAYER

SIROLIMUS – ACTIVE PHARMACEUTICAL INGREDIENT

Sirolimus is an immunosuppressive agent. Sirolimus inhibits T-lymphocyte activation and proliferation that occurs in response to antigenic and cytokine (Interleukin [IL]-2, IL-4, and IL-15) stimulation by a mechanism that is distinct from that of other immunosuppressants. In cells, sirolimus binds to the immunophilin, FK Binding Protein-12 (FKBP-12), to generate an immunosuppressive complex. This complex binds to and inhibits the activation of the mammalian Target Of Rapamycin (mTOR), a key regulatory kinase. This inhibition suppresses cytokine-driven T-cell proliferation, inhibiting the progression from the G1 to the S phase of the cell cycle. The sirolimus is intended to reduce restenosis as ancillary medicinal substance to coronary intervention using the Ultimaster Tansei Sirolimus Eluting Coronary Stent.

DEGRADABLE POLYMERS – INACTIVE INGREDIENTS

Ultimaster Tansei Sirolimus Eluting Coronary Stent has a coating consisting of two layers: a primer layer and a drug matrix layer. The primer layer and drug carrier polymer, poly (D,L-lactide-co-caprolactone) copolymer is expected to degrade within 3-4 months. The drug coating is applied abluminal, leaving the luminal side of the stent free from drug as such enhancing endothelial coverage.

BALLOON DELIVERY CATHETER

The delivery catheter is a high pressure, semi-compliant balloon delivery catheter that has two radiopaque markers, which fluoroscopically mark the ends of the stent to facilitate proper stent placement.

The active balloon length is closely sized to the length of the stent to prevent over-expansion of the tissue proximal or distal to the stent.

At the proximal end of the system is a female luer lock connector hub. This hub connects to the balloon inflation lumen. The guide wire enters the distal tip of the catheter and exits 25 cm proximal to the tip.

The surface of catheter is partially coated with hydrophilic polymer coating which generates lubricity when wet.

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ENGLISH

INDICATIONS

Ultimaster Tansei Sirolimus Eluting Coronary Stent System is indicated for improving myocardial blood flow in patients with stenotic lesions in coronary arteries, including but not limited to patients with STEMI, NSTEMI, acute coronary syndrome, diabetes mellitus, multivessel disease, bifurcation lesions, patients older than 65 years, male and female patients, patients with totally occluded lesions, long lesions, lesions residing in small coronary vessels, restenotic lesions including in-stent restenosis, ostial lesions, lesions in left main coronary artery. The Ultimaster Tansei stent system is suitable for both femoral and radial approach.

CONTRAINDICATIONS / RECOMMENDATIONS

Contraindications

- Patients in whom anti-platelet and/or anti-coagulant therapy is contraindicated.
- Patients with lesion(s) that prevents complete inflation of an angioplasty balloon.
- Patients with known allergy to L605 cobalt-chromium alloy.
- Patients with known allergy to nickel.
- Patients with known hypersensitivity to sirolimus or its structurally related compounds.
- Patients with known hypersensitivity to lactide polymers and caprolactone polymers.
- Patients with known hypersensitivity to contrast agent that cannot be controlled prophylactically prior to Ultimaster Tansei Sirolimus Eluting Coronary Stent implantation.
- Patients with extreme vessel tortuosity that may impair stent placement.

Recommendations

- It is strongly recommended not to implant stent in women who are pregnant.
- Effects of sirolimus during lactation have not been evaluated, therefore it is strongly recommended to avoid breast feeding when this stent is implanted.

WARNINGS / PRECAUTIONS

WARNINGS

- Judicious selection of patients is necessary since Percutaneous Coronary Intervention with the use of stents carries the risk of stent thrombosis, vascular complications and/or bleeding events. Hence patients should be maintained on clinically adequate post-procedural antiplatelet therapy (aspirin and thienopyridine, or appropriate antiplatelet agents).
- Only physicians who have received appropriate training should perform implantation of the stent.
- Any advancement after introduction of the delivery catheter into the vessel should be done under high resolution fluoroscopy. When resistance is felt during manipulation, determine the cause of the resistance before proceeding.
- Proper judgment is necessary to select lesion for direct stenting since insufficiently prepared lesion may lead to stent dislodgement.
- Ensure that the aluminium pack and blister pouch have not been damaged or opened as this may compromise the stability and the sterile barrier.

STORE THE DEVICE BETWEEN 1 – 30°C IN THE ALUMINIUM PACK.

The device is packed under oxygen free conditions.

Aluminium pack includes an oxygen absorber and a desiccant. Discard them without opening.

After opening the aluminium pack, use the device within 12 hours.

Do not store the device in the blister pouch.

PRECAUTIONS

Stent Handling – Precautions

- For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- Do not use a device that has reached or exceeded its expiry date.
- Ensure that the blister pouch has not been damaged or opened as this may compromise the sterile barrier.
- Use immediately after opening the blister pouch.
- The entire operation should be carried out aseptically.
- Do not use if the stent is exposed to abnormal rubbing or contact with objects other than the guiding catheter or opened hemostatic valve prior to implant.
- Do not rub or scrape the stent coating.
- Do not displace or remove stent on or from its delivery system as it may damage the stent and/or lead to stent embolization. Stent system is intended to perform as a system.
- Stent should not be used in conjunction with other delivery systems.
- Delivery system should not be used in conjunction with other stents.
- Special care must be taken not to handle or in any way disrupt the stent on the balloon. This is most important during catheter removal from packaging, catheter removal from holder, removal of protector sheath from stent, catheter placement over guide wire and advancement through rotating hemostatic valve adapter and guiding catheter hub.
- Do not “roll” the mounted stent with your fingers as this action may loosen the stent from the delivery balloon.
- Do not expose system to organic solvent. Use only the appropriate balloon inflation media. Do not use air or any gaseous medium to inflate the balloon as this may cause uneven expansion and difficulty in deployment of the stent.
- Do not attempt to straighten the proximal shaft (hypotube) as it may cause the catheter to break if it is accidentally bent.
- Exposing the stent to fluids before implantation is not recommended. Exposure to fluids prior to implantation may result in premature release of drug.

Stent Placement – Precautions

- Do not introduce negative pressure, or pre-inflate delivery system prior to stent deployment other than as directed.
- Always select an appropriate size of the stent as an undersized stent may result in inadequate expansion of the lesion while an oversized stent may lead to inadequate expansion of the stent or damage to the vessel wall.
- Always verify whether the stent is well apposed against vessel wall because incomplete stent apposition may lead to stent thrombosis.
- When treating multiple lesions in the same vessel, stent the distal lesion prior to stenting the proximal lesion. Stenting in this order avoids crossing the proximal stent with the distal stent and reduces the chances for dislodgement.
- Do not expand the stent if it is not properly positioned in the vessel. (See Stent System Removal – Precautions)
- Placement of a stent has the potential to compromise side branch patency.
- Do not exceed rated burst pressure as indicated on the device label. Use of pressure higher than specified may result in balloon rupture with possible intimal damage and dissection.
- An unexpanded stent may be retracted into the guiding catheter one time only. Subsequent movement in and out through the distal end of the guiding catheter should not be performed as the stent may be damaged or dislocated. In case of stent dislodgement, stent retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the coronary vasculature and/or the vascular access site. Complication may include bleeding, hematoma or pseudoaneurysm.

Stent / System Removal - Precautions

Stent introduction into the coronary artery is limited to one time only as dislodgement may occur.

Should unusual resistance be felt at any time during either lesion access or removal of the stent delivery system pre-stent implantation, carefully attempt to pull the stent delivery system back through the guiding catheter. If resistance is felt in doing so, or if resistance is felt during the removal of the stent delivery system post-stent deployment, the delivery system and guiding catheter must be removed as a single unit.

When removing the delivery system and guiding catheter as a single unit:

- Do not attempt to retract an unexpanded stent into the guiding catheter while engaged in the coronary arteries. Stent damage or dislodgement may occur.
- Position the proximal balloon marker just distal to the tip of the guiding catheter.
- Advance the guide wire into the coronary anatomy as far distally as safely possible.
- Tighten the rotating hemostatic valve to secure the delivery system to the guiding catheter, then remove the guiding catheter and delivery system as a single unit.
- Failure to follow these steps and/or applying excessive force to the delivery system can potentially result in stent dislodgement or damage to the stent and/or delivery system components.
- It is necessary to retain guide wire position for subsequent artery/lesion access, leave the guide wire in place and remove all other system components.

Post Implant - Precautions

- Care must be exercised when crossing a newly deployed stent with a coronary guide wire, IVUS catheter, OCT catheter, balloon or other stent delivery system to avoid disruption of the stent geometry.
- Patients should be maintained on clinically adequate post-procedural antiplatelet therapy (aspirin, thienopyridine or other appropriate antiplatelet agents) according to the current guidelines. In case of need, dual antiplatelet therapy can be discontinued earlier, but not before one month.
- Magnetic Resonance Imaging (MRI)
Non-clinical testing has demonstrated the Ultimaster Sirolimus Eluting Coronary Stent, overlapped configuration (max. 2 stents x max. OD4mm x 38mm = 73.6mm total length) is MR conditional.

It can be scanned safely under the following conditions:

MR Conditional for

- static magnetic field of 1.5 Tesla and 3 Tesla only, with
- spatial gradient field of 36 T/m and less
- spatial gradient field product of 99 T²/m and less
- theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of < 2 W/kg at 1.5 Tesla, (related to 3.6°C temperature increase; allowable level in accordance with CEM43 concept), 73.6 x 4.0 mm, overlapped configuration
- < 2 W/kg at 3 Tesla, (related to 3.6°C temperature increase; allowable level in accordance with CEM43 concept), 73.6 x 4.0 mm, overlapped configuration for 15 minutes of continuous MR scanning.

Temperatures and SAR have been derived from computer modeling with human realistic anatomy (no cooling effects considered).

In non-clinical worst-case phantom testing the Ultimaster Sirolimus Eluting Coronary Stent, overlapped configuration (max. 73.6 x 4.0 mm) produced a temperature rise of less than 5.2°C (with a background temperature increase of ~1.1°C) at a maximum whole body averaged specific absorption rate (SAR) of ~2.3 W/kg assessed by calorimetry for 15 min. of continuous MR scanning with whole body coil in a 1.5 Tesla Intera, Philips Medical Systems (Software: Release 12.6.1.4, 2012-11-05) MR Scanner.

In non-clinical worst-case phantom testing the Ultimaster Sirolimus Eluting Coronary Stent, overlapped configuration (max. 73.6 x 4.0 mm) produced a temperature rise of less than 10.1°C (with a background temperature increase of ~2.4°C) at a maximum whole body averaged specific absorption rate (SAR) of ~2.1 W/kg assessed by calorimetry for 15 min. of continuous MR scanning with whole body coil in a 3 Tesla Magnetom Trio, Siemens Medical Solutions (Software: Numaris/4, syngo MR A30) MR Scanner.

The Ultimaster Sirolimus Eluting Coronary Stents have not been tested in simultaneous combination with other devices.

MR image quality is compromised if the area of interest is in the same area or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this implant.

The results referenced above are obtained from Ultimaster Sirolimus Eluting Coronary Stent evaluation tests. The stent of Ultimaster Tansei Sirolimus Eluting Coronary Stent System is identical to the Ultimaster Sirolimus Eluting Coronary Stent.

INDIVIDUALIZATION OF TREATMENT

The risks and benefits of Sirolimus-eluting stent should be considered for each patient before (using) implanting Ultimaster Tansei Sirolimus Eluting Coronary Stent. Physicians are responsible for assessing patient appropriateness for undergoing stent implantation prior to procedure.

OPERATOR'S MANUAL

Inspection Prior to Use

- Carefully inspect the stent delivery system package for damage to the sterile barrier. Prior to using the Stent System, carefully remove the system from the package and inspect for bends, kinks, and other damage.

Materials Required / Quantity Material.

- Appropriate guiding catheter (minimum inner diameter of guiding catheter is 1.42mm (0.056"))
- 2 – 3 syringes (10 – 20 ml)
- 1,000 u/500 ml Heparinized Normal Saline (HepNS)
- 0.36 mm (0.014") x 175 cm (minimum length) guide wire
- Rotating hemostatic valve with appropriate minimum inner diameter [2.44 mm (0.096")]
- Diluted contrast medium 1:1 with heparinized normal saline (HepNS)
- Inflation device
- Pre-deployment dilatation catheter
- Three-way stopcock
- Torque device
- Guide wire introducer
- Appropriate arterial sheath
- Appropriate anticoagulation and antiplatelet drugs

Preparation

Guide wire Lumen Flush

Step Action

1. Carefully remove the stent system from its holder, then remove the stent protector sheath from over the stent.

CAUTION Carefully slide off the protector sheath from the stent by fixing the sheath at the distal end between the thumb and finger while gently pulling sheath and attached stylet.

2. Verify that the stent is centered on the balloon and located between the radiopaque balloon markers.

CAUTION Do not use if any defects are noted.

3. Flush the guide wire lumen with HepNS using the flushing needle supplied with Ultimaster Tansei stent system. Insert the flushing needle into the tip of the catheter and flush until solution exits the guide wire port.

CAUTION Avoid manipulating the stent while flushing the guide wire lumen as this may dislocate the stent on the balloon.

Delivery Procedure

Step Action

1. Prepare vascular access site according to standard practice.
2. Predilate lesion with PTCA catheter.
3. Remove the PTCA catheter.
4. Open rotating hemostatic valve on the guiding catheter as widely as possible.
5. Backload Delivery System onto proximal portion of guide wire while maintaining guide wire position across target lesion.

CAUTION Confirm that the guide wire OD does not exceed 0.36 mm (0.014"). If a large size guide wire has been used, exchange the wire in the standard manner.

6. Advance stent Delivery System over guide wire to target lesion. Utilize radiopaque balloon markers to position stent across lesion; perform angiography to confirm stent position.

CAUTION Take care not to damage the delivery catheter and stent when advancing the delivery catheter over the guide wire.

Deployment Procedure

Step Action

1. Before deployment reconfirm the correct position of the stent relative to the target lesion via the catheter markers.
2. Attach the inflation device to the delivery catheter hub and apply negative pressure to purge the balloon of air.
3. Under fluoroscopic visualization, inflate the balloon to at least nominal pressure and maintain for 15-30 seconds to deploy the stent but do not exceed the labeled rated burst pressure (see label on packaging or enclosed compliance sheet).
4. Optimal expansion requires the stent to be in full contact with the artery wall, and with the stent internal diameter matching the size of the reference vessel diameter.
5. Stent wall contact should be verified through routine angiography or intravascular ultrasound.
6. Deflate the balloon by pulling a vacuum with the inflation device. Make sure the balloon is fully deflated before any attempted movement of the catheter.
7. Confirm adequate stent expansion by angiographic injection through the guiding catheter.

Removal Procedure

Step Action

1. Ensure balloon is fully deflated.
2. Fully open rotating hemostatic valve.
3. While maintaining guide wire position, withdraw Stent Delivery System.
Note: Should unusual resistance be felt at any time during either lesion access or removal of the stent delivery system before stent implantation, the entire system should be removed. See stent/System Removal – Precautions section for specific stent delivery system removal instructions.
4. Tighten rotating hemostatic valve.
5. Repeat angiography to assess stented area. If an adequate expansion has not been obtained, exchange to a balloon catheter of appropriate balloon diameter to achieve proper stent apposition to the vessel wall. Do not dilate the stent beyond the limit tabulated below.

Nominal stent diameter (ID)	Post-Dilatation limits (ID)
Φ2.25mm, 2.5mm, 2.75mm, 3.0mm	4.5mm
Φ3.5mm, 4.0mm	5.5mm

6. Final stent diameter should match reference vessel diameter.

INSTRUCTION FOR SIMULTANEOUS USE OF TWO DEVICES IN GUIDING CATHETER (KISSING BALLOON TECHNIQUE)

6Fr Compatibility – Any combination of one Ultimaster Tansei Sirolimus Eluting Coronary Stent System (I.D. 2.25mm – 4.0mm) and one Hiryu® (2.25mm-4.00mm) or one Accuforce® (2.00mm-4.00mm) PTCA balloon catheter can be used simultaneously within a 6Fr (I.D.:1.8mm) guiding catheter.

The technique can be performed as per the instructions listed below:

1. Insert Ultimaster Tansei Sirolimus Eluting Coronary Stent System using the instructions provided.
2. Insert one Hiryu or one Accuforce balloon catheter, track to the target site and inflate the balloon.
3. Removing the catheters: Remove one catheter and its associated guide wire completely prior to removing the other catheter and its associated guide wire.

CAUTION Care should be taken when introducing, torquing and removing one or both devices to avoid entanglement.

DRUG INTERACTION

Drugs that act through the same binding protein (FKBP) may interfere with the efficacy of sirolimus.

Sirolimus is metabolized by CYP3A4. Strong inhibitors of CYP3A4 (e.g. ketoconazole) might cause increased sirolimus exposure to levels associated with systemic effects, especially if multiple stents are deployed. Systemic exposure of sirolimus should also be taken into consideration if the patient is treated concomitantly with systemic immunosuppressive therapy. Based on the results of the human pharmacokinetic study, the systemic effect of sirolimus after single stent implantation is considered negligible.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with coronary stent placement include but are not limited to:

- Abrupt vessel closure
- Acute myocardial infarction
- Allergic reaction to anti-coagulation and/or anti-thrombotic therapy, contrast material, or stent and/or delivery system materials or any other PCI mandatory medication
- Aneurysm
- Arrhythmias, including ventricular fibrillation and ventricular tachycardia
- Arteriovenous fistula
- Cardiac tamponade
- Cardiogenic shock
- Death
- Emboli, distal (air, tissue or thrombotic emboli)
- Emergent Coronary Artery Bypass Surgery
- Failure to deliver the stent to the intended site
- Hemorrhage, requiring transfusion
- Hypotension / Hypertension
- Infection and pain at insertion site
- Ischemia, myocardial
- Myocardial infarction
- Nausea and vomiting
- Prolonged angina
- Pseudoaneurysm
- Renal failure
- Restenosis of stented segment
- Rupture of native and bypass graft
- Stent compression
- Stent embolization
- Stent migration

- Fever
- Hematoma
- Thrombosis (acute, subacute, or late)
- Total occlusion of coronary artery
- Unstable or stable angina pectoris
- Stent thrombosis / occlusion
- Stroke / cerebrovascular accident
- Vessel dissection
- Vessel perforation
- Vessel spasm

Potential adverse events that may be associated with sirolimus drug and polymer coating. Sirolimus administration is limited to intra-coronary stent delivery. Therefore adverse events are not fully characterized but are considered to be consistent to those noted in sirolimus oral administration including:

- Abnormal liver function tests
- Anemia
- Arthralgias
- Changes in lipid metabolism which may include hypertriglyceridemia or hypercholesterolemia
- Diarrhea
- Hypersensitivity to the drug (sirolimus or its excipients) or to the polymer (or individual components) including anaphylactic / anaphylactoid type of reactions
- Hypokalemia
- Immune suppression, especially in patients with hepatic insufficiency or who are taking medications that inhibit CYP3A4 or P-glycoprotein
- Infections
- Interstitial lung disease
- Leukopenia
- Lymphoma and other malignancies
- Myalgia
- Thrombocytopenia

Because of the low systemic exposure to sirolimus after stent implantation, it is very unlikely that any of the adverse events (apart from hypersensitivity reaction) associated with oral administration of sirolimus, will occur.

HOW SUPPLIED

STERILE AND NON PYROGENIC in undamaged and unopened blister pouch. This device is sterilized by e-beam.

CONTENTS: One Terumo balloon expandable, Sirolimus Eluting coronary stent mounted on a rapid exchange delivery system. One Flushing needle.

STORAGE: Store between 1 and 30°C.

DISPOSAL: After use dispose of delivery system in accordance with local regulations.

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CONTRE-INDICATIONS / RECOMMANDATIONS

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All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners.
**Ashitaka Factory of Terumo Corporation**